



## FOLDING BARIATRIC PORTERING CHAIR: MR4588

MR



### Dimensions:

#### Overall:

Height: 925mm  
Width: 742mm  
Depth: 1055mm  
Width (folded): 475mm

#### Seat:

Height: 550mm  
Width: 630mm  
Depth: 550mm

This new chair has been designed specifically to transport bariatric patients to and from the MRI suite and is MR safe.

With the increase of obesity becoming more of an issue in recent years, this chair allows patients to be transported safely to an open scanner.

The MR4588 is approximately 30% wider than our standard portering chair and has been weight tested to 300kg (approx 47 stone) static weight.

With no magnetic parts included in the design, the chair can transport the patient in to the scanner room.

Armrests raise and lower to below seat level to ease patient transfer.

**THE MR4588 IS FULLY  
TESTED MR SAFE**

### Wardray Premise Ltd

Hampton Court Estate  
Summer Road  
Thames Ditton  
Surrey  
Tel: +44 (0) 20 8398 9911  
Fax: +44 (0) 20 8398 8032  
E-mail: sales@wardray-premise.com

Full details of all products and services can be viewed on our website:

[www.wardray-premise.com](http://www.wardray-premise.com)

QUALITY WITHOUT  
COMPROMISE SINCE 1909