

Magnetic Testing of Products at 7.0T

ITEM: MR2520 – Utility Trolley with Drawers

Manufactured by Wardray Premise Limited at:
Northern Regional Office
Unit 2 - Springvale Works
Elland Road
Brighouse
West Yorkshire
HD6 2RN
UK

This item was tested using the 7.0T Siemens Terra Scanner at the University of Cardiff – Brain Research Imaging Centre (CUBRIC) on 1st July 2019. The item was taken right up to the magnet bore and turned around to present different sides of the item to the magnet.

The item was not attracted to the magnet.

CE

Justine Colquhoun

.....
Justine Colquhoun
Director

Date: 22.4.2020