



MRI REVERSIBLE PRISM GLASSES

These reversible Prism Glasses are specifically designed for the MRI patient and has the advantage of being adjustable to either feet or head view direction.

A clear view of the scanner bore can considerably reduce anxiety and movement due to claustrophobia. Can be used with either plastic arms or head strap adjustment for extra comfort. Glasses can be stored in a box for hygiene and safety.

Features include:

- ▶ adjustable viewing direction up or down
- ▶ adjustable for adult and child patients
- ▶ can be used with MRI viewing devices such as fMRI screen or monitor
- ▶ MR safe to 7 Tesla



MR5650: MRI Reversible Prism Glasses



QUALITY WITHOUT COMPROMISE
FROM THE UK' S LEADING
RADIATION SHIELDING

WARDRAY PREMISE LTD

Hampton Court Estate
Summer Road
Thames Ditton
Surrey KT7 0SP
Tel: +44 (0) 20 8398 9911
Fax: : +44 (0) 20 8398 8032
E-mail: sales@wardray-premise.com

Full details of all products and services can be viewed on our website:

www.wardray-premise.com

SUCCESSFULLY TRADING
FOR 100 YEARS

GEN.S.11/230