

**Magnetic Testing of Products at 7.0T**

ITEM: MR1505 – Cube Trolley

Manufactured by Wardray Premise Limited at:  
Northern Regional Office  
Unit 2 - Springvale Works  
Elland Road  
Brighouse  
West Yorkshire  
HD6 2RN  
UK

This item was tested using the 7.0T Siemens Terra Scanner at the University of Cardiff – Brain Research Imaging Centre (CUBRIC) on 1<sup>st</sup> July 2019. The item was taken right up to the magnet bore and turned around to present different sides of the item to the magnet.

The item was not attracted to the magnet.

**CE**

*Justine Colquhoun*

.....  
Justine Colquhoun  
Director

Date: 22.4.2020