



WARDRAY  
PREMISE

# MR5502: Adjustable Height Patient Trolley

The MR5502 adjustable height patient trolley has been designed specifically for use in MR departments & is tested to work safely in MRI environments up to 7 Tesla.

The MR5502 has brand new & improved features to improve both the patient & users experience.

## 300kg Weight Capacity:

A robust 300Kg weight limit accommodates a broad range of patients, prioritising their comfort and safety

## Central & Directional locking castors:

All 4 castors are controlled by a simple foot operated brake bar for efficient central & directional locking





**WARDRAY  
PREMISE**

# MR5502: Adjustable Height Patient Trolley

## Cross Contamination:

Zero Velcro on mattress offers improved cleaning standards



## Improved Mobility:

Larger twin wheeled castors & adjustable push handles ensures maximum trolley manoeuvrability



## Trolley Dimensions

Length	2010mm
Width (cotsides down)	785mm
Width (cotsides up)	730mm
Maximum Height (cotsides down):	925mm
Minimum Height (cotsides down):	625mm
Maximum Load	300kg

## Trendelenburg:

Gas spring assisted Trendelenburg ensures smooth & controlled transitions

