



MR BEB'X PAEDIATRIC SECURING SYSTEM

Beb'x Paediatric Securing System: MRX EIM 106

This securing system is compatible with all medical imaging techniques including MRI, CT, X-ray and Ultrasound (e.g of the hip) all in the horizontal position.

This new system has been developed for young children up to 10kg. It is simple to set up and due to its soft mattress and securing strips very comfortable for the child. Thanks to the air mattress the patient has an optimal comfort position which makes the medical examination less stressful.

The patient securing systems are made by using the new super soft 'magic strip' straps. These straps have superseded the Velcro® method of securing in line with infection control procedures.

This system permits the child to feel safe, comfortable and totally relaxed.

- ▶ 1 x patient board tray
- ▶ 1 x inflatable mattress
- ▶ 4 x 'magic strip' straps
- ▶ 1 x pair of wrist straps



Beb'x Paediatric Securing System:
MRX EIM 106



MRX EIM 107

Beb'x Rotating Paediatric Securing System: MRX EIM 107

This unit has the same features as the 106 model but enables the child to be scanned with 3 positions of angle 0°, 45° and 90°.

There are no screws or buttons. A simple adjustment is required to set the height and any angle of rotation

- ▶ 1 x patient board tray with rotation system
- ▶ 1 x inflatable mattress
- ▶ 4 x 'magic strip' straps



Beb'x Rotating Paediatric Securing System:
MRX EIM 107

QUALITY WITHOUT COMPROMISE
FROM THE UK'S LEADING
RADIATION SHIELDING EXPERTS

WARDRAY PREMISE LTD

Hampton Court Estate
Summer Road
Thames Ditton
Surrey
KT7 0SP

Tel: +44 (0) 20 8398 9911
Fax: +44 (0) 20 8398 8032
E-mail: sales@wardray-premise.com

Full details of all products and services can be viewed on our website:

www.wardray-premise.com

More information overleaf

MR BEB'X PAEDIATRIC SECURING SYSTEM

General Information

The following information relates to both the MRX EIM 106 and MRX EIM 107 models.

Cleaning

For both units, we recommend using tap water, followed by the spraying of a water based solution and wiping with disposable alcohol free solution disinfectant wipes.

All cleaning must be followed by complete drying. Support bands and wrist straps are machine washable at a maximum temperature of 30°. Do not tumble dry.

Additional accessory options

MRX FMS 278	Magic strips on a roll of 20m x 60mm
MRX CIM 359	Packet of 5 inflatable mattresses
MRX FMS 301	Air pump for mattress inflation
MRX FMS 299	Pair of wrist straps

Technical Data

- ▶ safety weight limit is 10kg.
- ▶ magic strips are 2000mm x 60mm in length.
- ▶ patient tray is 840mm (L) x 340mm (W).
- ▶ weight of empty patient tray is 2kg.
- ▶ board material is polycarbonate 6mm thick.
- ▶ storage temperature: +5°C to +40°C.
- ▶ humidity: 10% to 90%

Classification: **CE 93/42/CEE: Class 1**



Beb'x Paediatric Securing System:
MRX EIM 106



Beb'x Rotating Paediatric Securing System:
MRX EIM 107

WARDRAY PREMISE LTD



SUCCESSFULLY
TRADING FOR
100 YEARS

Hampton Court Estate
Summer Road
Thames Ditton
Surrey
KT7 0SP

Tel: +44 (0) 20 8398 9911
Fax: +44 (0) 20 8398 8032

E-mail: sales@wardray-premise.com