



Mobile Screens Panoramic & Standard

A large selection of mobile screens are manufactured with panoramic (full width) windows or standard smaller windows. White or light grey plastic laminate finishes are available from stock. These screens are particularly useful in departments where space may be an issue.

A wide range of other colours are available to special order. The enamelled legs fitted with easy running anti-static castors provide excellent stability.

Panoramic Mobile Screens

Ref:	Lead	Overall		Window	
		Height	Width	Height	Width
PMS 800/150	1.5mm	2000mm	850mm	800mm	800mm
PMS 800/200	2.0mm	2000mm	850mm	800mm	800mm
PMS 1000/150	1.5mm	2000mm	1050mm	800mm	1000mm
PMS 1000/200	2.0mm	2000mm	1050mm	800mm	1000mm
PMS 1200/150	1.5mm	2000mm	1250mm	800mm	1200mm
PMS 1200/200	2.0mm	2000mm	1250mm	800mm	1200mm

The example sizes listed are 2000mm high overall with glass glazing 800mm high



Standard Mobile Screens

Ref:	Lead	Overall		Window	
		Height	Width	Height	Width
SMS 900/150	1.5mm	2000mm	900mm	400mm	400mm
SMS 900/200	2.0mm	2000mm	900mm	400mm	400mm
SMS 1000/150	1.5mm	2000mm	1000mm	400mm	400mm
SMS 1000/200	2.0mm	2000mm	1000mm	400mm	400mm

The example sizes listed are 2000mm high overall with windows 1500mm from floor to centre sight.